

COLD STERILANT®

SAFE COLD STERILISATION & DISINFECTING OF MEDICAL EQUIPMENT

YOU CAN CLEAN WITHOUT DISINFECTING, BUT YOU CANNOT DISINFECT WITHOUT CLEANING.

COLD STERILANT® containing Quadex® (V6) is a HIGH-LEVEL STERILANT AND DISINFECTANT, containing anti-rust technology, for the cold sterilising of surgical instruments and heat sensitive equipment. It is a Bactericidal, Fungicidal, Yeasticidal, Sporocidal when used undiluted (ready to use) on pre-cleaned surgical instruments. Use in industrial and domestic areas.

DIRECTIONS FOR USE:

After the pre-cleaning step (**INSTRUMENT PRE-CLEANER**), submerge pre-cleaned instruments/equipment into undiluted (ready to use) suspension and:

- **FOR INTERMEDIATE LEVEL DISINFECTING:**
Leave for 5 minutes.
- **FOR HIGH-LEVEL DISINFECTING/STERILISING:**
Leave for 30 minutes.

Quadex® (V6) exhibits the latest features desired of a modern disinfectant.

- NON-TOXIC
- NON-STAINING
- NON-IRRITANT
- NON-CORROSIVE
- pH NEUTRAL
- FAST-ACTING
- NON-HAZARDOUS
- BROAD-SPECTRUM
- EFFICACY
- WATER BASED
- REMOVES BIOFILM
- ENVIRONMENTALLY SOUND
- NEUTRALISES (organic and inorganic odours and malodours)
- RESIDUAL (30-day laboratory tested residual - unlike most disinfectants that merely kill and evaporate)

Cold Sterilant® is available in:

- 1 Litre - Bottle
- 5 Litre - Bulk

Ambit of use:

- B/S C/S



"ABSOLUTE PEACE OF MIND"

Cell: 086 111 2661 • sales@quadex.co.za • www.quadex.co.za

