

STERIFOG® BULK

HIGH-LEVEL DISINFECTING FOGGING SOLUTION (ULV)

YOU CAN CLEAN WITHOUT DISINFECTING, BUT YOU CANNOT DISINFECT WITHOUT CLEANING.

STERIFOG® containing Quadex® (V6) is a HIGH-LEVEL DISINFECTANT that treats exposed and hard to reach areas and is safe to use where people and patients congregate. Can be used in operating theatres, clinics, medical and dental surgeries, veterinary surgeries, mortuaries, microbiological laboratories, ICUs, burn units, patient rooms and the aesthetics industry. Odour-neutralising, air-freshening and anti-allergenic 7-in-1 formulation.

DIRECTIONS FOR USE:

Apply by using ULV fogging machine to disinfect hard surface areas. For heavily soiled areas a preliminary deep cleaning is required. Leave for 30 minutes to high-level disinfect. Area/work surfaces should be visibly clean.

Quadex® (V6) exhibits the latest features desired of a modern disinfectant.

- NON-TOXIC
- NON-STAINING
- NON-IRRITANT
- NON-CORROSIVE
- pH NEUTRAL
- FAST-ACTING
- NON-HAZARDOUS
- BROAD-SPECTRUM
- EFFICACY
- WATER BASED
- REMOVES BIOFILM
- ENVIRONMENTALLY SOUND
- NEUTRALISES (organic and inorganic odours and malodours)
- RESIDUAL (30-day laboratory tested residual - unlike most disinfectants that merely kill and evaporate)

Sterifog® is available in:

- 20 Litre - Bulk (RTU)

Fragrances:

- Plain and Vanilla

Ambit of use:

- C/S D/S

GMDN Code:

- 47631



"ABSOLUTE PEACE OF MIND"

Cell: 086 111 2661 • sales@quadex.co.za • www.quadex.co.za

