

STERIFOG® SPRAY

HIGH-LEVEL DISINFECTING SPRAY

YOU CAN CLEAN WITHOUT DISINFECTING, BUT YOU CANNOT DISINFECT WITHOUT CLEANING.

STERIFOG® SPRAY containing Quadex® (V6) is a HIGH-LEVEL DISINFECTANT, safe to use where people or animals congregate. Due to its low level of toxicity, it can be used in operating theatres, clinics, medical and dental surgeries, veterinary clinics, microbiological laboratories and patients' rooms. Odour-neutralising, air-freshening and anti-allergenic 7-in-1 formulation.

DIRECTIONS FOR USE:

Shake can well. Spray STERIFOG® into the air and on surfaces. STERIFOG® has been developed to create CLEAN AIR for the benefit of asthma sufferers and the removal of harmful irritants and toxins present in indoor atmospheres. Leave for 5 minutes to high-level disinfect. Area/work surfaces should be visibly clean.

Quadex® (V6) exhibits the latest features desired of a modern disinfectant.

- NON-TOXIC
- NON-STAINING
- NON-IRRITANT
- NON-CORROSIVE
- pH NEUTRAL
- FAST-ACTING
- NON-HAZARDOUS
- BROAD-SPECTRUM
- EFFICACY
- WATER BASED
- REMOVES BIOFILM
- ENVIRONMENTALLY SOUND
- NEUTRALISES (organic and inorganic odours and malodours)
- RESIDUAL (30-day laboratory tested residual - unlike most disinfectants that merely kill and evaporate)

Sterifog® Spray is available in:

- 350ml - Cans

Fragrances:

- Plain, Vanilla and Sunshine

Ambit of use:

- C/S D/S

GMDN Code:

- 47631



"ABSOLUTE PEACE OF MIND"

Cell: 086 111 2661 • sales@quadex.co.za • www.quadex.co.za

